



MULTI TREATMENT CHAIR

MULT-V2-60/70-150



(Picture of the model with options: pair of security side rails and reclining foot-rest)

DETAILS :

- Multi-treatment chair with electrical height adjustment from 56 to 96 cm
- Frame comprised of 2 telescopic columns providing the couch high stability
- Shaped base in very resistant plastic : Optimal hygiene, light grey color
- Hand remote with fixation on the right or left side of the chair
- One section mattress with a high comfort, width 600mm or 700mm (in option), made of visco-elastic foam, M1 class, 18 colors available for upholstery
- Back rest and leg rest adjustment via 2 independent electric jacks
- Inclination of the seat is electric
- Trendelenburg position of the chair
- Equipped with 2 side rails that can be used for several accessories and allow to adjust the arm rest lengthwise, equipped in standard with exam sheet roller, adaptable head rest and push bar in the back rest
- Operation: Equipped with 4 directional double running design wheels diam. 150mm, with a central brake system used with pedals on both wheels in the front and brake bar at the back of the chair
- Dimensions: Upholstery Length. 2000mm x Width 600 mm (or 700mm optional) x Height from 56 to 96 cm
Length of the back rest= 900 mm / Length of the seat rest= 540 mm/
Length of the leg rest= 540 mm
Width overall = 850mm

OPTIONS :

- **ACC-125** : 4 wheels diameter 125cm
- **NR** : Adjustable black foot rest, adjustable with gas jack
- **ACC-FIXE** : 4 adjustable fixed blades
- **ACC-CAPOTMULTIK** : Painted hood MULTIKA Version 2, following base range
- **ACC110POLY** : Polyurethane blood test splint with serrated block system but without clamps, the unit

- **ACCBATTERIE01** : Battery for chair MULTIKA
- **ACCOVALIA03**: Pedal OV2 – FMV2 – MULT-V2
- **ACC514-MULTIK02** : Plastic protection for leg rest
- **ACCBARRLITMUL02** : Adjustable security side rails
- **ACCCLAMEDIAM16**: Pair of clamps available on stainless steel accessory rail 25x10
- **ACCFEMINA02** : Pair of stirrups for MULTIKA without clamps
- **ACCFEMINA964** : Pair of leg rests fitted on kneecap available on polyurethane foam
- **SLFEMIREPOSE-UR** : Cover for leg rest ACCFEMINA964, the unit
- **ACCMULTIK02** : Emergency position button
- **ACCMULTIK05** : Reading tablet



PICTURES OF THE PRODUCT IN DIFFERENT POSITIONS :

