

A blue surgical chair is the central focus, positioned in an operating room. The chair has a high backrest, adjustable armrests, and a reclining table. It is mounted on a white base with four wheels. In the background, there is a whiteboard, a monitor, and other medical equipment. The floor is a light blue color.

NovyMed

Surgery chair

DOC Neo

The next generation

novymed.com

Mobilizing 
Healthcare

The next generation surgery chair: **DOC Neo**



The possibility to bring the chair in a **completely horizontal position**.



250 KG
safe working load

Side rails at leg-seat and back section.



Standard lithium battery (LFP)
Eight times longer lifespan than lead-acid versions and better for the environment.

Future proof controls with the possibility to increase the number of memory positions and CPR-function added.


Foldable arm rests enables free lateral movement of the patient in case of emergencies.

Push handles

Increased maximum height adjustment with a lowest position of 54 cm up to a maximum height of 86 cm.

The chair is constructed with **ready and available components**, easily checked during maintenance and repairs; parts efficiently replaced when necessary.

Completely redesigned driving mechanism for improved and unmatched manoeuvrability!



NovyMed Export BV
Vaart 7C, 4206 CC Gorinchem
The Netherlands

+31 (0)85 077 33 64
export@novymed.com
novymed.com

Mobilizing Healthcare